Two simple inoculation methods for Phytophthora blight and root rot studies

Mike Benson
NC State University
Raleigh

Steve Tjosvold
University of California
Cooperative Extension,
Watsonville

NC STATE Ornamental Pathology





Phytophthora inoculum

- Rice grain method
 - 25 g long grain rice
 - 18 ml deionized water
 - Autoclave 40 min

Holmes, KA and Benson, DM. 1994 Plant Disease 78:193-199.








Rice grain method

Transfer agar disks





How many rice grains?

























Foliage symptoms: root rot







Scale: 1 = Healthy; 4 = Dead

Root symptoms









Year to year: consistency

2006

| | | Foliar rating | | Top wt | Root | |
|----------------------------|-------------|--------------------|---------|---------|---------|------------------|
| Treatment and Rate/100 gal | Application | Jul 20 | Aug 18 | Sep 7 | (oz) | rot ^y |
| Untreated control | | 2.4 b ^z | 2.9 bc | 2.9 b | 0.59 d | 2.9 b |
| Aliette 80WP 80 oz | Spray | 1.1 c | 1.3 f | 1.4 e | 2.09 ab | 1.0 c |
| Non-inoculated control | | 1.0 c | 1.6 def | 1.6 cde | 1.74 ab | 1.0 c |

2004

| | Application | Disease severity ^z | | Top wt | Root |
|--|-------------|-------------------------------|--------|---------|---------------------------|
| Treatment and rate /100 gal applied at a 30-day interval | | 22 Jul | 09 Sep | OZ | rot (1-5) ^y |
| Untreated control | | 2.4 a × | 3.0 a | 0.7 e | 3.0 a |
| Aliette 80W 80.0 oz | Spray | 1.0 b | 1.2 c | 2.2 abc | 1.2 c |
| Non-inoculated control | | 1.2 b | 1.0 c | 2.3 abc | 1.0 c |

Hinodegiri azalea

Rice grain no. & Disease

Hinodegiri azalea: June 6 – September 19, 2006

| | Top wt. | Root ro | | |
|---------------|---------|---------|--|--|
| Treatment | (g) | (1-5) | | |
| Untreated | 29 a | 1.2 b | | |
| 1 rice gr/pot | 25 a | 2.6 ab | | |
| 3 rice gr/pot | 25 a | 1.8 ab | | |
| 6 rice gr/pot | 18 ab | 2.0 ab | | |
| 9 rice gr/pot | 7 b | 3.2 a | | |



Pot volume = 2.5 L

Uniform inoculum placement

- Rice grain particles
 - Pre- & post-emergence damping-off studies





Pulverize & Sieve






No. 10 sieve with 2-mm-openings



Weigh inoculum

- How much to use?
 - Trial & error







Seed & Cover








Pre-emergence damping-off



Did you notice the rice color change?







Foliar inoculation chamber

- Steve Tjosvold
 - University of California Cooperative Extension, Watsonville







Chamber close-up

- Inoculum
 - Mycelial disk
 - Spore suspension
 - 5,000 zoospores/cap





Chamber in place









Uses

- Pathogenicity tests
- Fungicide trials
- Environmental studies
- Hosts/cultivar evaluations



Questions?







