

Quantitative multiplex detection of (plant) pathogens, including *Phytophthora* species, based on PRI-lock probe technology and the OpenArray platform

Phytophthora, Pythium workshop ICPP, August 2008

Peter Bonants, Ronald van Doorn, Odette Mendes, Richard van Hoof and Cor Schoen



PLANT RESEARCH INTERNATIONAL
WAGENINGEN UR

Introduction

- Detection and identification in disease management strategies
 - Fast
 - Accurate
 - Sensitive
 - **Multiplex**

- Current technologies
 - Low level of multiplexing
 - Low throughput
 - Laborious
 - **Not always quantitative**



Introduction


- Many targets to be detected simultaneously

- Pathogens

- Fungi
- Oomycetes
- Bacteria
- Nematodes
- “Viruses”

- Beneficial microorganisms


Viruses.....



Fungi.....



Nematodes...



Bacteria...




Introduction


- For multiplex detection there is an increasing interest in quantification of the different organisms
- What are the possibilities for quantitative multiplex detection?
 - Sensitive DNA quantification only with qPCR
- Quantification in PCR for **multiple targets** is theoretically problematic



Padlock Probe Principle




Padlock Probe Principle




Multiplex detection in recirculation water

- *Phytophthora* sp
- *P. nicotianae*
- *A. tumefaciens*
- *V. dahliae*




Phytophthora micro-array


- PLP for *Phytophthora* sp
- PLPs for different *Phytophthora* species
- ligation with mixture of PLPs
- amplification with generic primers
- detection on *Phytophthora* micro-array



PRI-lock Probe Principle



PRI-lock




Hybridization
Ligation

Ligation Dependent



Quantitative Multiplex Target Detection

PRI-lock ligation followed by singleplex amplification




Quantitative Multiplex Target Detection

Multiplex PRI-lock ligation followed by **singleplex** amplification



PRI-lock ligation product



Pre-spotted primer combinations



Universal qPCR analyses



Quantitative Multiplex Target Detection

- Testing PRI-lock/primer specificity in a multiplex (96 well) system.


Target template	PRI-lock probe	<i>Myr. ror.</i> primers	All <i>Phyt.</i> primers	<i>Phyt. inf.</i> primers
<i>Myrothecium roridum</i>	<i>M. roridum</i>	21	40	40
<i>Phytophthora infestans</i>	All <i>Phytophthora</i>	40	21	40
<i>P. Infestans</i>	<i>P. infestans</i>	40	40	22
<i>M. roridum</i>	PRI -lock mixture	21	40	40
<i>P. infestans</i>	PRI -lock mixture	40	21	22
<i>M. ror. + P. inf.</i>	PRI -lock mixture	21	21	22
No target	PRI -lock mixture	40	40	40

- PRI-lock/primer combinations are specific
- Ct values of the PRI-locks are not influenced by the presence of other templates and PRI-locks in the mixture



Quantitative Multiplex Target Detection


Real-time quantification in a Biotrove 'PCR array'




Quantitative Multiplex Target Detection

- 48 subarrays on each plate
- Each subarray consists of 64, 33 nL through-holes in an 8x8 pattern
- Primer pairs spotted in the through-holes
- 3072 reactions in each array
- Loading by capillary action

Through-hole cross-section








Quantitative Multiplex Target Detection





Loading three OpenArray plates to an OpenArray Cycler







Quantitative Multiplex Target Detection

Real-time quantification in a Biotrove 'PCR array'



Primer pair	All Phyt	P. inf	M. ror	M. hapla	C. spat	E. car	R. sol 2.2	V. albo/ tri	V. dahl	R. sol 4.1	G. Proteo
Targets											
P. inf/ M. ror	17.2	17.5	18.1	ND	ND	ND	ND	ND	ND	ND	ND
M. hapla/ C. spa/ E. car	ND	ND	ND	19.2	19.4	17.4	ND	ND	ND	ND	19.7
V. dahl/ Rhiz 4-1/ Rhiz 2-2/ V. tri	24.6	ND	ND	ND	ND	ND	17.3	19.4	17.6	19.3	ND
P. inf/ M. hapla/ M. ror/ V. dahl/ V. tri/ E. car	17.0	17.3	18.5	19.3	ND	17.6	ND	19.2	18.1	ND	20.4
MO	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND



Quantitative Multiplex Target Detection

- Testing PRI-lock/primer specificity in the Biotrove OpenArray platform

Primer pair	Phyt spp.	<i>P. inf</i>	<i>M. ror.</i>	<i>M. hap.</i>	<i>A.tum.</i>	<i>E. car.</i>	<i>F. oxy.</i>	<i>V. alb./ tri.</i>	<i>V. dah.</i>	<i>R. sol. 4-1</i>	<i>R. sol. 4-2</i>
Targets											
<i>P. infestans</i>	16.4	15.4	--	--	--	--	--	--	--	--	--
<i>A.tumefaciens</i>	--	--	--	--	16.2	--	--	--	--	--	--
<i>A. tum. / F. oxy. / V. dah.</i>	--	--	--	--	16.2	--	17.7	--	15.9	--	--
<i>P. inf. / M. ror. / M. hap. / A. tum. / E. car. / F. oxy. / V. alb. / V. dah. / R. sol. 4-1 / R. sol. 4-2</i>	15.9	15.0	17.0	16.9	15.9	15.4	17.3	18.0	15.8	14.4	18.8
Milli Q water	--	--	--	--	--	--	--	--	--	--	--

- PRI-lock/primer combinations are specific
- Ct values of the PRI-locks are not influenced by the presence of other templates in the mixture



Quantitative Multiplex Target Detection

■ Advantages:

- Quantitative
- High specificity
- Single and multiplex target detection independent
- Target recognition and amplification independent
- Universal TaqMan / SYBR Green PCR conditions
- High-throughput
- Low background

■ Disadvantages:

- Low copy numbers of ligated PRI-locks in nanoliter wells



Conclusions

- Working PRI-lock probe multiplex detection system:
 - Currently 30 PRI-lock probes
- Fields of application:
 - Multiplex quantitative target (pathogen) detection
 - Microbial community analysis
 - Multiplex quantitative gene expression analysis
- Instrumentation:
 - Standard real-time PCR machines
 - Nanoliter 'PCR arrays' – OpenArray platform (BioTrove)



Acknowledgements

- Plant Research International B.V.
Wageningen, The Netherlands
 - Marianna Szemes
 - Carolien Zijlstra
 - Richard van Hoof
 - Ronald van Doorn
 - Els Nijhuis
 - Odette Mendez
 - Cor Schoen

- Isogen Life Science
IJsselstein, The Netherlands
 - Nicole Wellens

- BioTrove Inc.
Woburn, USA
 - Jonathan Schimmel
 - Elen Ortenberg

References

- Szemes et al. Nucleic Acids Research, 2005, Vol. 33, No. 8 e70
 - Van Doorn et al. BMC Genomics, 2007, 8:276
-



Thank You For Your Attention

© Wageningen UR



PLANT RESEARCH INTERNATIONAL
WAGENINGEN UR